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
The screenshot shows the ODFPO1eFigure2 software interface. At the top, three ODF plots are displayed for different reflections: $\{0,0,2\}$ (1211.6), $\{1,0,2\}$ (922.6), and $\{1,0,3\}$ (1632.8). The main control panel includes the following sections:

- Files select:** ASC(RINT-PC) and 002-ref.ASC 102-ref.ASC 103-ref.ASC
- Calculation Condition:** Previous/Next buttons, C:\CTR\DATA\Mg#002-ref.ASC, hkl: 0,0,2
- Background delete mode:** DoubleMode (checked), SingleMode, LowMode, HighMode, Nothing. BG defocus: DSH1.2mm+Schulz+RSH5mm.
- Smoothing:** + α (checked), 3, Arithmetic mean, Disp.
- Peak slit:** 7.0 mm, BG Slit: 7.0 mm, PeakSlit / BGS: checked, BG Scope: 80.0 deg, 90.0 deg.
- AbsCalc:** Ref, Trans, Schulz reflection method, Absorption coefficient: 133.0 1/cm, Thickness: 0.2 cm, 2Theta: 34.54 deg, 1/Kt (checked), Profile.
- Defocus file Select:** Transmission defocus HKL+T (checked), Defocus(1) functions file, Defocus(3) function files folder (Calc unbackdefocus): BB106mm, Defocus(2) function files folder (Calc backdefocus): DSH1.2mm+Schulz+RSH5mm (circled in red), Limit Alfa Defocus value: Free(LimitValue=0.0), Select crystal Cubic (unchecked), 1/Ra (checked), Profile.
- Smoothing for ADC:** Cycles: 5, Weight: 8, Disp.
- Normalization:** CTR (checked), Connect, Average, OutFiles: Asc, M TexAsc, Ras, TXT, TXT2 (checked).
- Buttons:** Cancel, Calc, Connect, Exit&ODF, ODF, ValueODFVF-B, ValueODFVF-A.
- Status:** Select crystal : Hexagonal, 19/09/30.

A close-up view of the software control buttons. The 'Calc' button is circled in red.

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

A close-up view of the software control buttons. The 'ODF' button is circled in red.

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

Recaldefocusfile
C:\CTR\DATA\Mg\LaboTex\CW\caldefocus\DEFOCUS_F.TXT

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