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
Search thickness で厚さを計算

Absorption coefficient 8.69 1/cm Thickness 0.003 cm

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



The screenshot shows the PoleOrientationTools 1.13ST[19/10/31] by CTR interface. On the left is a circular pole figure plot with a color scale from blue to red. The main window contains several panels:



- File Help**: A list of operations with corresponding buttons:
 - TPF,HPF,TXT,TXT2 Full Polefigure TXT → Orientation → fnd,frd,ftd Display
 - TXT2 Reflection(001) PoleFigure → NDOrientation → fnd Display
 - TXT2 Reflection PoleFigure → CreateExpPoleFigure → fnd,frd,ftd Display
 - TXT2 PoleFigure → RandomLevelCalc → Cr Display
 - TXT2 Trans.Reflection polefigureTXT2 → PFConnection → Complete polefigure
- Files select**: ASC(RINT-PC) and 020_ref.ASC.
- Calcratation Condition**: Includes background delete mode (DoubleMode, SingleMode, LowMode, HighMode, Nothing), BG defocus (DSH1.2mm+Schulz+RSH5mm), Trans blinds angle (30.0), Peak slit (10.0 mm), BG Slit (10.0 mm), PeakSlit / BGS (80.0 deg, 90.0 deg), AbsCalc (Ref, Trans), Schulz reflection method, Absorption coefficient (0.01 1/cm), Thickness (0.1 cm), 2Theta (36.0 deg), and 1/Kt.
- Defocus file Select Transmission defocus HKL+T**: Options for Defocus(1) functions file, Defocus(3) function files folder (BB185mm), and Defocus(2) function files folder (LINE-BB-185mm).
- Smoothing for ADC**: Cycles (5), Weight (8), Disp.
- Normalization**: CTR, Connect, Average, Asc, MTexAsc, Ras, TXT, TXT2.
- Buttons**: Cancel, Calc, Connect, Exit&ODF, ODF, ValueODFVF..., ValueODFVF-A.
- Bottom status**: Select crystal : Orthorhombic 19/09/16.

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The screenshot shows the AbsCoefSearch 1.00ST[19/10/31] by CTR interface. On the left is a circular pole figure plot. The main window contains several panels:

- Material**: Material Polyethylene.TXT.
- Select TXT2**: TXT2 C:\VCTR\DATA\PE\020_ref_chB60_2.TXT.
- AbsCalc**: Schulz reflection method, Absorption coefficient (8.69 1/cm), Thickness (0.0098 cm), 2Theta (36.42 deg), 1/Kt. The thickness value 0.0098 is circled in red.
- Buttons**: TD, RD, Search thickness, AbscalcPoleFigure.

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