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The screenshot shows the ODFPolefigure1.5.1.63 software interface. The window title is "ODFPolefigure1.5.1.63 by CTR PDuser ctr helpertex". The interface includes a menu bar (File, Linear(3D), Toolkit, Help, InitSet, Rp%, Minimum, All background, Transmissionblinds=30.0) and a toolbar. The main area contains several panels: "Files select" (with a file browser icon), "Calculation Condition" (Previous, Next, hkl: 0,0,0), "Background delete mode" (DoubleMod, SingleMode, LowMode, HighMode, Nothing, Minimum($\alpha\beta$), MiniAver X, 0.5), "AbsCalc" (Ref, Trans, Schulz reflection method, Absorption coefficient: 133.0, Thickness: 0.2, cm, 2theta: 0.0, deg, 1/Kt), "Defocus file Select" (Normalization, Polynomialdegr.: 0, TenckhoffFitting, TXT2), "Smoothing(for ADC)" (Cycles: 2, Weight: 15, Afterconnection), "Normalization" (CTR, Connect, Average), "CenterData", "Search minimum EqualAngleRp%(Cubic only)", "OutFiles" (Ras, Asc, TXT2, TXT), "Cancel", "Calc", "Connect", "ODF File", "Select crystal: Monoclinic 19/12/02").


Annotations with arrows point to specific features:

- Red box around "Files select" and "Background delete mode" with arrow pointing to "データ入力" (Data input).
- Red box around "Background delete mode" with arrow pointing to "バックグラウンド" (Background).
- Green box around "AbsCalc" with arrow pointing to "defocus".
- Yellow box around "Smoothing(for ADC)" with arrow pointing to "平滑化" (Smoothing).
- Blue box around "Normalization" with arrow pointing to "規格化" (Normalization).
- Yellow box around "Smoothing" and "RD" with arrow pointing to "平滑化" (Smoothing).
- Yellow box around "RD" with arrow pointing to "RD 補正" (RD correction).
- Yellow box around "OutFiles" with arrow pointing to "Rp%改善" (Rp% improvement).

The screenshot shows the ODFPoleFigure2 3.94 software interface. The window title is "ODFPoleFigure2 3.94 by CTR PDuser ctr helpertex". The interface includes a menu bar (File, Linear(absolute)3D, Toolkit, Help, InitSet, BGMode, Measure, Condition Free, OverlapRevision, MinimumMode, Rp%, Normalization) and a toolbar. The main area contains several panels: "Files select" (with a file browser icon), "Calculation Condition" (Previous, Next, hkl: 0,0,0), "Background delete mode" (DoubleMode, SingleMode, LowMode, HighMode, Nothing, Minimum($\alpha\beta$), MinimumAverage(α)X, 0.5, BG defocus: DSH1.2mm+Schulz+RSH5mm, Minimum mo...), "Peak slit" (7.0 mm), "BG Slit" (7.0 mm), "PeakSlit / BGS..." (checked), "BG Scope" (80 deg, 90 deg), "Set", "Disp", " α Inhibit", "Smoothing" (+ α : 3, Arithmetic mean, Disp), "RD" (0.0, Interpolation, Full, Disp), "AbsCalc" (Ref, Trans, Schulz reflection method, Absorption coefficient: 300.0, Thickness: 0.1, cm, 2Theta: 0.0, deg, 1/Kt), "Defocus file Select" (Transmission defocus, HKL+T, Defocus(1) functions file, Make defocus function files by TXT2, Files, Normalization, degree of a polynomial: 0, TenckhoffFitting, TextDisp, Defocus(3) function files folder(Calc unbackdefocus): BB185mm, Limit Alfa Defocus value: Free(LimitValue=0.0), Defocus(2) function files folder(Calc backdefocus): DSH1.2mm+Schulz+RSH5mm, Search minimum Equal Angle Rp%(Cubic only), 1/Ra, Profile), "Smoothing for ADC" (Cycles: 5, Weight: 7, Disp, After connection), "Normalization" (CTR, Connect, Average), "CenterData", "OutFiles" (Asc, MTexAsc, Ras, TXT, TXT2), "Cancel", "Calc", "Connect", "Exit&ODF", "ODF", "ValueODFVF-B", "ValueODFVF-A", "Select crystal: Monoclinic 19/12/02").

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